# RED DE OBSERVATORIOS MICROBIANOS

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APHA (American Public Health Association) (1998). Standard methods for the examination of water and wastewater. En: Clescerl LS, Greenberg AE, Eaton AD (eds) Book 20. APHA, Washington, DC, p 4-487.

Fernández Zenoff V., Siñeriz F. & Farías M.E. (2006). Diverse responses to UV-B radiation and repair mechanisms of bacteria isolated from high-altitude aquatic environments. *Applied and Environmental Microbiology*, 72(12), 7857-7863. <https://doi.org/10.1128/AEM.01333-06>

Lozada M., Zabala M.S., García P.E., Diéguez M.C., Bigatti G., Fermani P., Unrein F., Dionisi H.M. (2022). Microbial assemblages associated with the invasive kelp *Undaria pinnatifida* in Patagonian coastal waters: structure and alginolytic potential. *Sci. Total Environ.* 154629. <https://doi.org/10.1016/j.scitotenv.2022.154629>

Strickland J.D.H., Parsons T.R. (1972). A practical handbook of sea-water analysis, 2nd edn. J. Fish. Res. Board Can. 167, 311.

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---

## 4.

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- APHA (American Public Health Association) (1998). *Standard methods for the examination of water and wastewater*. En: Clescerl LS, Greenberg AE, Eaton AD (eds) Book 20. APHA, Washington, DC, p 4-487.
- Fernández Zenoff V., Siñeriz F. & Farías M.E. (2006). Diverse responses to UV-B radiation and repair mechanisms of bacteria isolated from high-altitude aquatic environments. *Applied and Environmental Microbiology*, 72(12), 7857-7863. <https://doi.org/10.1128/AEM.01333-06>
- Lozada M., Zabala M.S., García P.E., Diéguez M.C., Bigatti G., Fermani P., Unrein F., Dionisi H.M. (2022). Microbial assemblages associated with the invasive kelp *Undaria pinnatifida* in Patagonian coastal waters: structure and alginolytic potential. *Sci. Total Environ.* 154629. <https://doi.org/10.1016/j.scitotenv.2022.154629>
- Strickland J.D.H., Parsons T.R. (1972). *A practical handbook of sea-water analysis*, 2nd edn. *J. Fish. Res. Board Can.* 167, 311.

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# MICROBIAL OBSERVATORY NETWORK

Protocols for standardization and water sampling in Latin America

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**Cover Photo: Dr. Cristian D. Torres**

# **Microbial Observatory Network**

*Protocols for standardization and water sampling in Latin America*

*Fermani Paulina, Gereá Marina, Graziano Martín, Mateus-Barros Erick, Sabio y  
García Carmen, Sánchez María Laura, Schiaffino Romina*

*Compilers*

*Gerea Marina and Bernal María Carolina*

*Translators*

*2023*

*First Edition*

## **Summary**

Anthropic impacts and climatic factors act jointly on ecosystems as forcing agents of change, often leading to their deterioration. The effect of these agents depends on geographic location, climate and vegetation, among other natural and anthropogenic factors that characterize the studied environment. In order to understand the impact of these factors over aquatic and terrestrial ecosystems it is important to have continuous, standardized and extended information over time.

Within aquatic ecosystems, microbial assemblages play a fundamental role in the cycles of matter and energy, and also have the potential to be robust indicators of the ecological and sanitary status of the ecosystem, due to their sensitivity and rapid response to environmental changes. In this sense, long term monitoring of microbial assemblages could provide valuable information on the characteristics of water bodies and the changes they undergo.

Considering these premises, in the framework of the first meeting of the collaborative network on microbial aquatic ecology in Latin America ( $\mu$ SudAqua), held in La Paloma (Rocha, Uruguay) in December 2017, the creation of a set of observatories was decided, whose main objective is to evaluate how anthropogenic impact and climatic factors, in a latitudinal gradient, influence the structure and dynamics of the microbial community at the continental level in the diversity of aquatic environments in the region. Thus, the Latin American Microbial Observatory Network was born.

Microbial observatories are valuable tools that have been established around the world in order to have access to systemized information on microscopic communities. However, they are scarce in Latin America. In this Observatory Network, observatory sites corresponding to different aquatic environments (freshwater, marine, lentic, lotic, with different types of anthropic impact) were selected based on their accessibility in order to achieve continuous sampling.

Likewise, the observatory sites chosen contemplate a minimum bimonthly (preferably monthly) and simultaneous frequency, in the surface layer of the water column (euphotic), and the measurement of simple parameters, using protocols agreed upon by its members, which facilitate long-term continuity and the possibility of making comparative analyses between the different sites.

This **book** aims to list the established and standardized protocols for sampling physicochemical parameters and biological communities in water, and subsequent analyses in the laboratory, taking into account the different water bodies.

The main objective of this network is to strengthen the collaborative links between research laboratories from different Latin American countries in order to enhance the capabilities and knowledge of each group. In this way, we intend to promote the exchange of knowledge to find collective answers to questions that concern us as a region.

## ***Preface***

***Fermani Paulina, Gereá Marina, Graziano Martín, Sabio y García Carmen, Schiaffino Romina***

Microbial observatories are a valuable tool to obtain systematized and organized information on microbial communities. Despite this, there are few observatories in Latin America that contemplate long-term studies, mainly due to the difficulty of maintaining the necessary infrastructure for this type of activity. In this sense, during the first meeting of the Aquatic Microbial Ecology Network ( $\mu$ SudAqua, <https://microsudaqua.netlify.app/>) held in December 2017 in Rocha (Uruguay), the need to create a Latin American Aquatic Microbial Observatory Network emerged. There, the long-term objectives for the Network were discussed and agreed upon. Thus, it was established as a general objective to form a set of observatories (sites) to evaluate how anthropogenic impact and climatic factors in a latitudinal gradient affect the structure and dynamics of the microbial community at the continental level in the diversity of aquatic environments in the region. We also hoped that this would allow us to generate a collective work among Latin American groups to strengthen regional and collaborative studies.

For the initial formation of the Latin American Aquatic Microbial Observatory Network, at the first meeting of the  $\mu$ SudAqua Network, eight observatory sites of different types of aquatic environments (freshwater, marine, lentic, lotic) throughout Latin America were proposed. In order to ensure the continuity of the sampling, the selection of sites strongly pondered the accessibility of the sites and that they should preferably be part of local projects that would help sustain the collection of samples. These sites were then re-evaluated based on the difficulties encountered during the Network's consolidation phase. Thus, the micro-observatory samplings contemplated a monthly frequency (possibly bimonthly, depending on site accessibility), in the surface layer of the water column (euphotic zone), and the measurement of simple parameters that facilitate long-term continuity. Subsequently, a survey of the proposed sites and the eventual incorporation of new observatories were carried out.

In order to address a global analysis of the parameters at each site, the common variables to be analyzed were discussed and defined during the first meeting in 2017. Once these parameters were established, a long process of protocol search for each parameter began which, within possibilities, could be used by all the observatories (rivers, streams, lakes, estuaries, etc.). Through surveys and virtual meetings, the equipment, techniques and protocols used by each site were evaluated and a consensus was reached on unified protocols for each of the proposed parameters. During this process, it became evident that the selected aquatic systems were very heterogeneous and the protocols that were useful for one type of system might not be useful for another system. For example, the Chlorophyll-a extraction efficiency with ethanol or acetone varies among environments (freshwater or marine). For this reason, it was defined to include both methods in this book. On the other hand, in the particular case of DNA extraction, sites that were already analyzing environmental DNA prior to the establishment of the network used different procedures (e.g. commercial kit or homemade technique; extraction with SDS and CTAB, etc.) and no common criterion was found to choose a single protocol, so it was decided not to establish a particular methodology. However, the present book includes one of the most widely used methodologies in current observatory sites for DNA extraction. All this information was compiled in a single document and shared with the observatory groups in order to initiate the Network's pilot sampling (April, 2019), with standardized protocols which allowed to reveal and discuss their advantages and difficulties during the initial period, in the second meeting of the observatory network. In November 2019, the second meeting of the  $\mu$ SudAqua was carried out in Chascomús (Buenos Aires, Argentina), where the second workshop of the microbial observatory network took place. There we shared the preliminary results of the first months of sampling of each observatory, and we moved forward in the discussion of different aspects of the network, sharing improvements and difficulties and establishing new goals. Until then, the Network included observatories located in Costa Rica, Brazil, Uruguay and Argentina, whose microbial communities and main limnological characteristics were monitored bimonthly. The observatory network included at that time eight different aquatic ecosystems (period 2019-2020): 4 lakes/lagoons, 2 rivers, 1 estuary and 1 marine site. Based on the discussion and due to

the impossibility of unifying the protocols among different aquatic sites, the main objective was redefined, and we decided to maintain a book of protocols that included only the most common techniques among the observatories. In this document, procedures for most of the parameters are unified, while in some cases more than one methodological option is present, depending on the type of environment analyzed (freshwater or marine) and others are left to the discretion of each group (e.g. environmental DNA). In addition, considering methodological uncertainties of some protocols, explanatory videos of some of the techniques were recorded and included in the document.

In 2020, as a consequence of the Covid-19 pandemic, most of the observatory sites interrupted the samplings for at least one and a half years. After a virtual meeting of the network in October 2021, we decided to re-start the synchronized sampling in April 2022, and we encouraged the incorporation of new observatory sites in the following stage. At present, the network includes 13 active observatory sites consisting of 3 lakes/lagoons, 4 rivers, 2 estuaries and 4 marine sites (<https://microsudaqua.netlify.app/>).

We presented the book of standardized protocols of the Observatory Network as a guide agreed upon at the 2019 meeting in Chascomús, which is currently used by all the observatory sites for sampling.

2023, December

## ***Microbial Observatory Network Participants (2017-2023)***

### **Argentina**

Allen Dohle, Sharon - INIBIOMA (UNComahue-CONICET)

Baliña, Sofía - IEGEBA (CONICET-UBA)

Barrena, Maité - CIMAS-CONICET

Bastidas Navarro, Marcela - INIBIOMA (UNComahue-CONICET)

Bernal, María Carolina - IEGEBA (CONICET-UBA)

Burgueño, Giuliana; CIMAS-CONICET, ESCiMar, UNCo

Cetra, Nicolás; CIMAS-CONICET, ESCiMar, UNCo

Fermani, Paulina - IBIOMAR-CONICET

García, Patricia - INIBIOMA (UNComahue-CONICET)

Gerea, Marina - INIBIOMA (UNComahue-CONICET)

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Protocol of Cytometry: Fernando Unrein and Andrea Malitz

Protocols of Dissolved Nutrients: Martín Graziano, Leonardo Lagomarsino

Protocols of Total Nutrients: María Luz Padulles and Eddy Gómez

Protocol of Environmental DNA: Carmen Sabio y García, Luciana Griffero, Paulina Fermani

Protocol of Dissolved Organic Matter and Dissolved Organic Carbon: Marina Gereá

## **Table of contents**

<b>General indications for sampling</b> .....	1
<b>Section 1: Preparation of material for sampling</b> .....	2
1.1 Chlorophyll-a and nutrients.....	3
1.2 Quantitative Phytoplankton .....	3
1.3 Cytometry.....	4
1.4 Environmental DNA.....	5
1.5 Dissolved Organic Carbon (DOC) and Dissolved Organic Matter (DOM).....	6
<b>Section 2: Field Material and sampling protocol</b> .....	7
2.1 Field measurements.....	8
2.2 Chlorophyll-a and nutrients.....	8
2.3 Quantitative Phytoplankton .....	8
2.4 Cytometry .....	8
2.5 Environmental DNA.....	10
2.6 Dissolved Organic Carbon (DOC) and Dissolved Organic Matter (DOM).....	10
<b>Section 3: Laboratory analysis protocols</b> .....	11
3.1 Turbidity.....	12
3.2 Chlorophyll-a and Nutrients.....	12
3.2.1 Determinación de Clorofila-a.....	13
3.2.2 Nutrientes disueltos.....	15
3.3 Quantitative Phytoplankton.....	15
3.4 Cytometry.....	16
3.5 Environmental DNA.....	17
3.6 Dissolved Organic Carbon (DOC) and Dissolved Organic Matter (DOM).....	20
<b>Section 4: References</b> .....	23
<b>Section 5: Epilogue</b> .....	24

---

## General indications for sampling

This section indicates how sampling will be carried out, how often, the sampling localization in the water body and what parameters will be measured.

### Sampling:

- ✓ Preferentially monthly sampling (eventually bimonthly).
- ✓ Morning sampling, during the second half of the month. If bimonthly, during even months (February, April, June, August, October, December).
- ✓ Sampling in the photic layer of the water body.
- ✓ How to identify the samples: TWO LETTERS FOR SITE\_yymmdd

### Parameters to measure:

#### 1) Chlorophyll-a (Chl-a)

#### 2) Nutrients:

- Total Nutrients ---> Total Phosphorus (TF), Total Nitrogen (TN)
- Dissolved Nutrients ---> Ammonium ( $\text{NH}_4^+$ ), nitrates ( $\text{NO}_3^-$ ), phosphates ( $\text{PO}_4^+$ )

#### 3) Quantitative Phytoplankton

#### 4) Cytometry

- Heterotrophic Bacteria
- Autotrophic picoplankton and nanoplankton

#### 5) Environmental DNA

#### 6) Dissolved Organic Matter (DOM) and Dissolved Organic Carbon (DOC)

## Preparation of material for sampling

The materials to consider for sampling are summarized in the following list. In the following sections the particularities of each parameter are mentioned:

- ✓ 250 or 500 or 1000 ml Plastics flasks (to determine total nutrients, dissolved nutrients/Chlorophyll-a, phytoplankton, environmental DNA, DOC/DOM). For Chlorophyll-a and environmental DNA in marine sites, collect samples in 5-liter drums.
- ✓ 10% HCl
- ✓ Distilled or ultrapure water (ej. Milli-Q®)
- ✓ Acidified Lugol solution (1%) (phytoplankton fixative)
- ✓ 5 or 4.5 ml sterile cryovials with specific fixative defined for each water body (to analyze prokaryotes, picoeukaryotes and nanoplankton by flow cytometry)
- ✓ Cryovial boxes (cytometry)
- ✓ Liquid N<sub>2</sub> thermos (cytometry)
- ✓ Pipette and tips
- ✓ Niskin or Van Dorn Bottle (depending on site)
- ✓ 50  $\mu$ m mesh

Particularities for each parameter will be detail in the following paragraphs.

## 1.1 Chlorophyll-a and Nutrients

Chlorophyll-a and Dissolved Nutrients: wash the plastic flasks with 10% HCl , then rinse with distilled water or ultrapure water at least 3 times. The size of the flask will depend on the volume to be filtered for the determination of Chlorophyll-a concentration. The filtered water (filtrate) will be used to determine the concentration of dissolved nutrients.

Total Nutrients: wash the 250 or 500 ml plastic flasks with 10% HCl , then rinse with distilled water or ultrapure water at least 3 times.

## 1.2 Quantitative Phytoplankton

Wash the 250 or 500 ml plastic flasks with distilled water or ultrapure water and prepare the acidified Lugol solution.

### Preparation of Acidified Lugol solution:

- ✓ Prepare 200 ml of Lugol 13,6% (30 grams of solute in 220 ml of solution = 13,6%*m/v*): 20 g IK + 10 g I<sub>2</sub> in 200 ml of distilled water (we recommend to use a magnetic shaker to accelerate the dissolution) + 10% of glacial acetic acid (20 ml).

Note 1: The acidified Lugol solution is stored in the dark and in the refrigerator (4°C), wrapping the bottle with aluminum foil.

Note 2: Lugol solution is also available in a ready-made format, which must be acidified with 10% glacial acetic acid.

### 1.3 Cytometry

The materials needed to prepare the samples for cytometry are listed below, depending on the fixative to be used (P+G, Glutaraldehyde, or GlyTE) and the type of organisms to be preserved (bacteria, picoautotrophs and nanoplankton):

- Label 4.5 ml cryovials and add 400  $\mu$ l of the fixative you will use (bacteria, picoautotrophs and nanoplankton).

Store the stock of cryovials with the aliquoted fixative in a  $-20^{\circ}\text{C}$  freezer.

The fixative preparations are listed below, according to the type of organisms to be preserved and the sampling site.

#### Preparation of P+G:

Prepare under a fume hood.

To prepare 1 liter of P+G (final concentration 1% Paraformaldehyde + 0.05% Glutaraldehyde):

- Heat 800 ml of PBS 1X to  $60^{\circ}\text{C}$  and maintain at that temperature with a magnetic shaker with heating.
- Add 100 g of Paraformaldehyde.
- Add NaOH lentils one at a time and continue to stir until the solution clarifies (should take 1-5 min).
- Remove from heat
- Add 20 ml of 25% Glutaraldehyde.
- Adjust to pH = 7.2 with HCl.
- Bring to 1 liter with 1X PBS.
- Filter through 0.22  $\mu\text{m}$  polycarbonate filter.
- Aliquot 400  $\mu$ l into each cryovial.
- Store in the freezer at  $-20^{\circ}\text{C}$  (preferably at  $-80^{\circ}\text{C}$ ).

Note: once thawed, the fixative should be stored at  $4^{\circ}\text{C}$  and used within one week, otherwise discarded!

Preparation of GlyTE (solution of glycerol + buffer TE):

- Prepare 100 ml EDTA solution 0.5 M:
  - Add 18.6 g EDTA + 80 ml ultrapure water (While mixing with the magnetic shaker, adjust to pH = 8 by slowly adding NaOH).
  - Add ultrapure water until a final volume of 100 ml.
  - Filter through 0.22  $\mu$ m polycarbonate membrane.
  
- Prepare 100X TE (1 M Tris; 100 mM EDTA):
  - Add 12.1 g of TRIS + 20 ml of 0.5 M EDTA (previously prepared)
  - Bring to 80 ml with ultrapure water (while mixing with the magnetic shaker, adjust to pH = 8 by slowly adding 37% HCl).
  - Add ultrapure water until a final volume of 100 ml.
  - Filter through 0.22  $\mu$ m polycarbonate membrane.
  
- Preparation of GlyTE stock:
  - Mix:
    - 20 ml 100x TE pH = 8,0
    - 60 ml of ultrapure water
    - 100 ml of molecular grade glycerol (use syringe)
  - Filter the GlyTE stock through a 0.22  $\mu$ m membrane and store in refrigerator (4°C).
  - Freshly prepared and filtered GlyTE is aliquoted into cryovials. The volume that remains unaliquoted should be re-filtered before use.

Note: autoclave the remaining 0.5M EDTA and 100X TE solutions and store until the next GlyTE preparation.

#### 1.4 Environmental DNA

Wash 250, 500 ml or higher volume (if needed) plastic flasks with 10% HCl once. Then rinse with distilled water 2-3 times, and a final rinse with ultrapure water. The size of the bottle will depend on the volume that is subsequently filtered.

Wash a filtration system in the same manner as the bottles, for laboratory use.

*Note: in general, 100 to 300 ml of sample is filtered in freshwater systems. In hypereutrophic lagoons such as Chascomús, 30 to 50 ml are filtered. In the sea, 1 to 4 liters are filtered.*

### **1.5 Dissolved Organic Carbon (DOC) and Dissolved Organic Matter (DOM)**

Wash a 250 mL plastic flask once with 10% HCl, and rinse with distilled water 3 times before bringing to the field to collect the water samples.

Wash a glass filtration system with 10% HCl once, and rinse with distilled water 3 times and with ultrapure water once, then muffle at 440°C for at least 1 hour. On the other hand, follow the same washing protocol on 100 ml caramel-colored bottles (which will be used to store the sample for DOC and DOM determinations) and muffle them at 440°C for at least 1 hour. The filtration system and bottle will be used in the laboratory after sampling.

---

## 2.

### ***Material to bring to the field and sampling protocol***

The materials needed for the field are listed below:

- ✓ Field sheet
- ✓ pH meter
- ✓ Conductimeter
- ✓ Oximeter and temperature sensor
- ✓ Secchi disk
- ✓ Niskin or Van Dorn bottle (depending on site)
- ✓ Pipette and tips
- ✓ Marker pen
- ✓ 50  $\mu\text{m}$  mesh or similar
- ✓ Refrigerator, ice and/or ice pack
- ✓ Liquid N<sub>2</sub> thermos (for cytometry)
- ✓ 250 or 500 ml pre-washed plastic bottles (total nutrients, dissolved nutrients/Chlorophyll-a, phytoplankton, environmental DNA, DOC/DOM)
- ✓ Acidified Lugol solution and plastic Pasteur pipette (phytoplankton)
- ✓ 4.5 ml Cryovials with P+G (cytometry): bacteria, picoautotrophs and nanoplankton
- ✓ 4.5 ml Cryovials with GlyTE (cytometry): bacteria, picoautotrophs and nanoplankton
- ✓ Echo sounder or graduated rod for water level measurement
- ✓ Distilled water to wash field sensors (Conductimeter, pH meter, Oximeter) and paper to dry.

Procedure:

**2.1 Field measurements**

Register water temperature, pH, conductivity and dissolved oxygen using field parameter sensors and following the equipment indications.

Identify and survey relevant variables for each environment. For example, width, depth, current velocity, flow, type of surrounding vegetation, etc., can be measured in streams. In lakes and sea, transparency with Secchi disc, surrounding vegetation, depth, etc.

**2.2 Chlorophyll-a and Nutrients**

Collect unfiltered samples of water in the previously rinsed plastic bottles and transport the samples in cold (4°C) and dark conditions.

**2.3 Quantitative phytoplankton**

Collect between 200 and 500 ml (depending on the study site) of unfiltered subsuperficial water in a previously rinsed bottle/flask. Fix sample with 1% acidified Lugol solution (final concentration). Keep in cold (4° C) and dark conditions until posterior analysis.

*Note: if quantification is not done right away, examine samples and add a drop of acidified Lugol (1%) so they do not lose coloration, and continue preserving.*

**2.4 Cytometry**

Collect water samples and pre-filter through a zooplankton net (50 µm mesh).

Distribute the sample into different cryovials containing specific fixatives according to the study site and microorganisms to be analyzed:

- Bacteria, picoautotrophs and nanoplankton: Add 4 ml of sample into 4.5 ml cryovials containing 400  $\mu$ l of either GlyTE or P+G.
- Homogenize each sample by inversion and keep in darkness for 10 minutes.
- Then, submerge samples in liquid nitrogen until they freeze and store in a box at  $-80^{\circ}\text{C}$  (otherwise store at  $-20^{\circ}\text{C}$ ).

Note: video-tutorial for processing and storage of samples available here:

[https://drive.google.com/file/d/1h6cbLzvIX\\_OJ5IX1of0xe03nkOW6sWsh/view?usp=sharing](https://drive.google.com/file/d/1h6cbLzvIX_OJ5IX1of0xe03nkOW6sWsh/view?usp=sharing)

Notes to consider for correct preservation of samples for cytometry:

- Fix the sample as soon as possible, avoiding the fixative to remain defrosted and at room temperature for long periods of time, particularly if using P+G.
- Ideally, take the liquid  $\text{N}_2$  thermos to the sampling site and fix the sample *in situ*.
- If not possible, collect the sample, carry it to the laboratory in cold ( $4^{\circ}\text{C}$ ) and dark conditions and fix it as soon as possible.
- Once defrosted, P+G must be used immediately to fix the sample. If not used right away, it can be stored at  $4^{\circ}\text{C}$  for 1 week at the most. After 1 week it should be discarded. It cannot be re-frozen nor kept at room temperature for long periods of time since the efficiency of fixation is altered.
- **Freezing** is a key point in the preservation of the sample. It is crucial that the fixed sample is frozen as soon as possible, and thus the necessity to count with

liquid N<sub>2</sub>. If you do not have liquid N<sub>2</sub> thermos, freeze at -80°C and, if this option is not available either, freeze at -20°C, although the latter is not recommended.

- In conclusion: **the lower the temperature => the faster the sample freezes => the better the quality of the sample!**
- **The sample must be kept frozen until analysis.** If transporting fixed samples, they must be frozen!

## 2.5 Environmental DNA

Rinse the previously washed bottles with water from the sampling site twice and collect the sample. It is recommended, especially in the more disturbed sites, pre-filtering through a 50 µm mesh.

Transport the sample to the laboratory in cold (4°C) conditions.

## 2.6 Dissolved Organic Carbon (DOC) and dissolved organic matter (DOM)

Collect 250 ml of water sample in the previously rinsed bottle (it is enough for both determinations). Transport the sample in dark conditions and ambient or cold temperature.

---

### 3.

#### ***Protocol for laboratory analysis***

The following list indicates the materials required for lab analysis posterior to sample collection:

- 10% HCl
- Distilled and/or ultrapure water
- Filtration system (Dissolved nutrients/Chlorophyll-a and DNA)
- Reagents, e.g. Hach (Dissolved nutrients)
- Potassium peroxodisulfate ( $K_2O_8S_2$ ) (Total nutrients)
- Boric Acid ( $H_3BO_3$ ) (Total nutrients)
- Sodium Hydroxide (NaOH) 1M (Total nutrients)
- 0.7-0.8  $\mu m$  pore glass fiber membranes (e.g. Whatman GF/F) muffled at 450°C for at least 4 hours (Dissolved nutrients/Chlorophyll-a)
- Aluminum foil (Chlorophyll-a)
- Glass filtering system washed with acid and muffled at 440°C for at least 1 hour (DOC and DOM).
- 0.22  $\mu m$  pore polycarbonate membranes (e.g. Millipore GTTP04700) (DNA)
- 0.3  $\mu m$  pore membranes(optional environmental DNA)
- Polyvinylidene fluoride filters (PVDF) or 0.22  $\mu m$  GV membranes (Millipore - GVWP04700) (DOC and DOM)
- 100 ml caramel-coloured glass bottle muffled at 440°C for at least 1 hour (DOC and DOM)
- Turbidimeter (if available)
- UV spectrophotometer with absorption spectrum capability

Procedure:

### 3.1 Turbidity

- In case you have a turbidimeter, save the specific volume of sample for this measurement according to the equipment requirements.

Note: In case you lack a turbidimeter, carry out a total suspended solids determination based on APHA protocols.

### 3.2 Chlorophyll-a and Nutrients

Total nutrients:

- Bottles containing samples collected in the field are preserved at -20°C until analysis.
- Dissolve 5 g of peroxodisulfate and 3 g of Boric Acid in 35ml of NaOH and bring to final volume of 100 ml with ultrapure water (peroxodisulfate is hard to dissolve; use a magnetic shaker).
- Once dissolved, use the solution in a proportion of 4 ml reagent : 30ml sample.
- Steam sterilize for 90 minutes.
- Cool and filter.
- Carry out the Soluble Reactive Phosphorus and Nitrate determinations.
- Determinations should be based on APHA or EPA protocols. For sea water, Strickland-Parsons (1972) methods are recommended.

Chlorophyll-a and dissolved nutrients:

- Immediately after the field sampling, filter a defined volume of sample through 0.7  $\mu\text{m}$  pore size muffled filters (e.g. Whatman GF/F). It is important to keep the material away from the light so as to reduce photooxidation. Filter sample until filter saturation. Filtered volume will depend on suspended materials. Thus, small volumes about 50-150 ml will be used in eutrophic environments while in oligotrophic or marine environments volumes can range from 1000 to 2000 ml.

- Registered the filtered volume.
- Place filter (folded in half, with retained material pointing inwards) in an aluminum foil envelope labeled with date, name of site, volume filtered and number of filters.
- Keep in  $-20^{\circ}\text{C}$  freezer to induce cell wall breakage and pigment release. Filters can be kept in these conditions for up to 3 weeks.

Nota 1: Filter will be used for Chlorophyll-a determination (Section 3.2.1).

Nota 2: Filtrate will be used for Dissolved nutrients determination (Section 3.2.2).

### **3.2.1 Chlorophyll-a determination:**

This step can be done a few days after filtering the sample and preserving the filter at  $-20^{\circ}\text{C}$ , but it is recommended that it is done before 3 weeks of storage have passed.

- Cut filters in small pieces (5 or 6 pieces) and place in flasks or tubes wrapped in aluminum foil. Add 8 ml of hot (between  $60$  and  $70^{\circ}\text{C}$ ) absolute ethanol per tube, as extraction solvent. In case 2 filters were used for a sample, add 12 ml of solvent instead of 8 ml.
- Vortex or shake vigorously and make sure all pieces of filter are submerged in the ethanol.
- Registered the volume of ethanol added. Absolute ethanol must be heated in an adequate container in a water bath. Control temperature with a thermometer placed in the ethanol.
- Leave samples resting for 24 hours in dark and cold ( $4^{\circ}\text{C}$ ) conditions to enable the extraction of photosynthetic pigments.

The day after the extraction:

- Centrifuge samples for 10 minutes at 3500-4000 rpm.

- Measure absorbance at 665 and 750 nm in the spectrophotometer. Previously measured pure ethanol as a blank.
- In the same cuvette, add 1 drop of HCl 1 N and, after a minute, repeat absorbance measurement at both wavelengths. Do not surpass this amount of HCl, since some accessory pigments may change their absorbance to that of pheophytin-a and thus interfere with the determination.

*Formula for the calculation of Chlorophyll-a concentration:*

$$[\text{Chlorophyll-a without phaeopigments}] = F [(Abs_{1_{665}} - Abs_{1_{750}}) - (Abs_{2_{665}} - Abs_{2_{750}})] k v$$

where Chlorophyll-a without phaeopigments is expressed in  $\mu\text{g}$  per liter; Abs1 = Absorbance before acidifying; Abs2 = Absorbance after acidifying; F = correlation factor that equates the reduction in absorbance with the initial Chlorophyll-a absorbance (2.43 for ethanol, 2.72 for methanol and 2.43 for acetone); k = coefficient of specific absorption (11.2 for ethanol, 11.62 for methanol and 10.48 for acetone); v = volume of extract in ml / (volume of sample x cuvette thickness in cm).

Considerations:

- ✓ Use ACS grade absolute ethanol.
  - ✓ Use cuvettes specific for chlorophyll. Wash with EXTRAN detergent and then rinse again with ethanol.
  - ✓ Calibrate with pure ethanol at 665 nm and 750 nm for the blank.
  - ✓ Repeat calibration every 4 samples to corroborate calibration.
  - ✓ Wash cuvettes with ethanol in between samples.
- Extraction with acetone for marine/saline water:
- Follow the same steps previously detailed, using acetone (non-heated) instead of hot ethanol as extraction solvent.

*Considerations:*

- ✓ Use 90% Acetone.
- ✓ Use cuvettes specific for chlorophyll. Wash with EXTRAN detergent and then rinse again with acetone.
- ✓ Calibrate with pure acetone at 665 nm and 750 nm for the blank.
- ✓ Repeat calibration every 4 samples to corroborate calibration.
- ✓ Wash cuvettes with acetone in between samples.

➤ Extraction with methanol for freshwater samples:

Alternatively to the use of hot ethanol, methanol can be used as an extraction solvent. In this case, a previous inter-calibration must be done in order to check the extraction capability of both solvents. Steps to be followed are the same as with pure ethanol.

**3.2.2 Dissolved nutrients:**

It is recommended that this determination is done the same day or immediate following days after the sample is collected. If nutrient measurement is not done the same day, the filtrate must be stored at -20°C until analysis.

*Note: parameter determination must be done according to APHA or EPA methods. We recommend: Ammonium, salicylate method (EPA Method 350.1) - alternative phenate method, APHA 4500-NH<sub>3</sub>-F; Nitrates: method of reduction by cadmium, APHA4500-NO<sub>3</sub>-E; Orthophosphate: method of vanadate-molybdate, APHA 4500-P-C - ascorbic acid method, APHA 4500-P- E. For seawater, we recommend following Strickland-Parsons (1972) methods.*

**3.3 Quantitative phytoplankton**

- Keep sample refrigerated (4°C) and in darkness until quantification using inverted microscope.

Note: Fixation with Lugol requires a periodic control of the coloration of stored samples, so as to replenish the Lugol which was oxidized. If samples appear discolored, add a drop of acidified Lugol.

### 3.4 Cytometry

By means of a flow cytometer, extracted and fixed samples for analysis of heterotrophic bacteria (1), autotrophic picoplankton (2) and nanophytoplankton (3) will be processed.

Due to the variability in flow cytometers available for each research group, it is not possible to standardize the settings. However, some criteria can be established for the analysis. In general terms:

- The *threshold* must always be set at the "FL1 or FITC" fluorescence for SybrGreen staining and at "FL3 or PerCP" for phytoplankton.
- Add 1  $\mu\text{m}$  sized Yellow-Green beads as a reference for size and fluorescence.
- Registered the analyzed and dilution volumes, flow velocity, tension, setting values and name of files.

#### 3.4.1 Heterotrophic bacteria:

- Bacteria are stained using SybrGreen I.

Preparation of working solution (WS) of SybrGreen I:

- > The 100X WS (working solution) of SybrGreen I is prepared diluting the commercial stock (10.000X) 1 : 100 in DMSO.
- For analysis in the flow cytometer, 100  $\mu\text{l}$  of sample are stained with 1  $\mu\text{l}$  of WS (1X final concentration).
- Incubate at room temperature for 10 minutes in darkness before analysis in the cytometer.
- Bacteria are typically analyzed in two biplots: FL1-SSC and FL3-FL1.

#### 3.4.2 Autotrophic picoplankton:

- Samples are analyzed in the cytometer without staining.
- Typically, FL3-SSC, FL3-FL2 and FL3-FL4 biplots are analyzed. The latter depend on the picocyanobacteria being rich in phycoerythrin (FL2) or phycocyanin (FL4).

#### 3.4.3 Autotrophic nanophytoplankton:

- The same sample and biplots as for autotrophic picoplankton are used, but settings are modified (voltage is lowered) so as to place picoeukaryotes in the lower-left section of the plot and better visualize the nanophytoplankton populations.

### **3.5 Environmental DNA**

#### ***Filtration and preservation procedure***

- Sterilize work surfaces and materials with 70% ethanol
- Filter between 100 and 300 ml of sample (except seawater: between 1 and 4 liters) through 0.22  $\mu\text{m}$  pore filters. To do this, a filtration system previously rinsed (once with 10% HCl, twice with distilled water and once with ultrapure water) will be used. Note: the filtrate obtained (between 100 and 150 ml) may be kept in a caramel-coloured bottle for posterior DOC analysis.
- Preserve each filter in a cryovial or sterile eppendorf tube without buffer, at  $-80^{\circ}\text{C}$ .

Note 1: *In severely disturbed sites a previous filtration through 3  $\mu\text{m}$  filters is recommended to avoid the posterior saturation of the 0.22  $\mu\text{m}$  filter with suspended particles.*

Note 2: *video-tutorial about filtration procedure available here:*

<https://drive.google.com/file/d/1h7XTJY8vcSX59xdISlInn6uTrGJpcbAK/view?usp=sharing>

### 3.5.1 Environmental DNA extraction procedure

We lack a standardized procedure, so each group uses their own protocols. For orientation purposes, we put at disposal two references of extraction procedures, one of them by means of a "home-made" method which is used by several groups (a) and the other using an extraction kit (b):

- a- Protocol of extraction with CTAB buffer: *Fernández Zenoff et al. (2006)*
- b- Protocol with extraction kit: *Lozada et al. (2022)*.

#### a- Protocol CTAB (Modified from *Fernández Zenoff et al. (2006)*)

##### **CTAB followed by extraction with solution of chloroform and isoamyl alcohol:**

- 1) Prepare CTAB (Cetyl Trimethyl Ammonium Bromide) lysis buffer solution (2% CTAB; 1.4M NaCl, 100 mM Tris-Cl pH 8, 20 mM EDTA pH 8). This solution is stored at room temperature and remains stable for several years. Before using it, add B-Mercaptoetanol at 0.2% final concentration and heat up to 60°C.
- 2) Aliquot CTAB solution in 2 ml eppendorf tubes and add to them each filter.
- 3) Incubate at 60 °C for 30 minutes.
- 4) Then, two purification steps are carried out: add 0.7 ml of chloroform-isoamyl alcohol (24:1) and centrifuge at 14,000 rpm for 10 minutes.  

Note 2: during centrifugation, make sure all tubes are oriented in the same direction, since pellets that will form are not visible with the naked eye.
- 5) Recover the aqueous phase in a 1.5 ml eppendorf tube (avoid catching pellet particles during this step).
- 6) Then, DNA is precipitated in cold isopropanol.

Note 3: Place in the fridge for at least 15 minutes. The longer it stays refrigerated, the more precipitation will occur. Consider leaving it in the fridge overnight.

7) Centrifuge at 14,000 rpm for 10 minutes

Note 4: during centrifugation, make sure all tubes are oriented in the same direction, since pellets that will form are not visible with the naked eye.

8) Remove the supernatant, carefully avoiding the pellet.

9) Then, rinsing steps are necessary. To do this, add cold 80% ethanol.

10) Mix and centrifuge for 2 minutes at 5,000 rpm.

Note 5: for oligotrophic environments, increase centrifugation time to 7 minutes.

11) Repeat rinsing with cold 80% ethanol.

12) Remove the supernatant, avoiding the pellet.

13) Drying and storage

Let the pellet dry at room temperature for 15 minutes with the opening of tubes pointing downwards so as to protect the pellet.

Note 6: drying can be done in a SpeedVac

14) Resuspend pellets in 100 $\mu$ L of ultrapure water.

Note 7: for oligotrophic environments, re-suspend in 50 $\mu$ L;

Nota 8: when adding the water, make sure it washes over the surface where the pellet is attached.

16) Determine DNA concentration in a Nanodrop®, Picodrop® or Fluorometer.

17) Dilute the DNA to a concentration of 10 ng/ $\mu$ L using ultrapure water.

18) Store in the freezer at -20°C.

### **b- Using extraction kits**

DNA can be extracted with the Fast DNA Spin Kit (MP Biomedicals Inc. USA) following the manufacturer's instructions, and quantification with fluorometer or Nanodrop®. For marine samples, procedures by Lozada et al. (2022) can be carried out.

### 3.6 Dissolved Organic Carbon (DOC) and Dissolved Organic Matter (DOM)

- Filter at least 200 ml of sample through 0.22  $\mu\text{m}$  PVDF or GV filters using a previously washed and muffled glass filtration system.
- Save the filtrate in two caramel-coloured bottles with teflon lids, previously washed and muffled. One of them will be destined to DOC analysis and the other to DOM analysis. Do not place aluminum foil between the bottle mouth and the lid.

*Note 1: use 0.22  $\mu\text{m}$  pore polyvinylidene fluoride filters (PVDF), GV Membrane filters (Millipore - GVWPO4700) or PES membrane filters, used previously for the DNA extraction.*

#### 3.6.1 Dissolved Organic Carbon (DOC):

- Preserve samples in cold (4°C) and darkness conditions.
- Process the sample in a carbon analyzer which has NPOC (non-purgable organic carbon) determination.

*Note 2: if the sampled site has a high load of organic matter (e.g. DOC  $\geq$  10 mg/L), add 3 drops of pure HCl for preservation.*

#### 3.6.2 Dissolved organic matter (DOM):

If possible, measure the absorption spectrum of the sample immediately after filtering. To do this, make sure you homogenize the temperature of the sample and of the ultrapure water which will be used as a blank in the spectrophotometer.

- Submerge the flasks containing the samples in a water bath at 20°C for 30 minutes before measuring. It is important to homogenize temperatures so as to avoid noise or irregular variations of the spectrum between 700 and 800 nm, a range that will then be used to correct the sample.
- Once the samples have reached 20°C, turn on the spectrophotometer and let its lamps warm up for at least 30 minutes.
- For the measurement, use a quartz cuvette with the appropriate light path length according to the type of environment (e.g. 1 cm, 5 cm or 10 cm).

- Configure the spectrophotometer to scan the spectrum of absorption of the sample between 200 and 800 nm. Extract the absorbance values for posterior use.
- **In order to optimize and homogenize calculations, the GESAP group (INIBIOMA-Bariloche) prepared an Excel form with macros for all members of the observatories network:**

<https://docs.google.com/spreadsheets/d/1nzTMZoUEo0AXijHO5dEfODJC1OMnaoyW/edit#gid=1978573616>

- This form automatically calculates the steps which are detailed next:
- Correct the absorption spectrum of the sample *a posteriori* of subtracting the absorption spectrum of the ultrapure water. To do this, average the absorbance values between 700 and 800 nm and subtract this average from each of the absorbances of the complete spectrum.
- Depending on your spectrophotometer, absorbance units (AU) should be converted into absorption coefficients (Abs) as follows:  $a = 2.303 A/L$ , where  $a$  is the Neperian absorption coefficient ( $m^{-1}$ ),  $A$  is the absorbance and  $L$  is the light path length through the cuvette (expressed in meters).
- Afterwards, using the corrected spectrum and the conversions into absorption coefficients, calculate:
  - a-  $a_{254}$
  - b-  $a_{350}$
  - c-  $a_{440}$
  - d- The spectral slope  $S_{275-295}$
  - e- SUVA 254 ( $a_{254} \cdot DOC$ )

f- SUVA 350 ( $a_{350}:\text{DOC}$ )

g-  $S_R$

- For each sample, record in the form provided by GESAP:
  - a- geographic location of sampling site,
  - b- volume of sample collected (mL),
  - c- filter used (PVDF, GV or fiberglass),
  - d- spectrophotometer used,
  - e- type of cuvette used (i.e., specify the light path length in meters, e.g. 0.01 m or 0.1 m).

## References

- APHA (American Public Health Association) (1998). Standard methods for the examination of water and wastewater. In: Clescerl LS, Greenberg AE, Eaton AD (eds) Book 20. APHA, Washington, DC, p 4-487.
- Fernández Zenoff V., Siñeriz F. & Farías M.E. (2006). Diverse responses to UV-B radiation and repair mechanisms of bacteria isolated from high-altitude aquatic environments. *Applied and Environmental Microbiology*, 72(12), 7857-7863. <https://doi.org/10.1128/AEM.01333-06>
- Lozada M., Zabala M.S., García P.E., Diéguez M.C., Bigatti G., Fermani P., Unrein F., Dionisi H.M. (2022). Microbial assemblages associated with the invasive kelp *Undaria pinnatifida* in Patagonian coastal waters: structure and alginolytic potential. *Sci. Total Environ.* 154629. <https://doi.org/10.1016/j.scitotenv.2022.154629>
- Strickland J.D.H., Parsons T.R. (1972). A practical handbook of sea-water analysis, 2nd edn. J. Fish. Res. Board Can. 167, 311.

---

## 5.

### *Epilogue*

Teamwork and collaboration are essential to carry out this great project. Undoubtedly, the contribution of each one of us and the predisposition to virtual and face-to-face meetings are the driving force to carry out with enthusiasm the temporary samplings and to plan the long-term projection of the Microbial Observatory Network. Here, we present some photos of the meetings and members together with some Latin American observatory sites.

*2017 - Teamwork of Microbial Observatory network (Uruguay)*


2019 - Teamwork of Microbial Observatory network (Argentina)

2021 - Teamwork of Microbial Observatory network (virtual meeting)

*2022 - Teamwork of Microbial Observatory network (Brasil)*

*Representative picture of marine observatory site, Punta Este, Puerto Madryn, Chubut, Argentina (Photo by Paulina Fermani)*

*Representative picture of lotic freshwater system observatory-site, Arroyo San Francisco, Claypole, Prov. Buenos Aires, Argentina (photo by Martín Graziano).*

*Representative picture of lentic freshwater system observatory-site, Laguna El Trébol, Bariloche, Prov. Río Negro, Argentina (Photo by Cristian D. Torres)*

*Representative picture of lentic freshwater system observatory-site, Reservorio Broa, Itirapina, San Pablo, Brasil (Photo by Erik Mateus-Barros)*